

M. Schumann<sup>1</sup>, R. Engels<sup>2</sup>, G. Kemmerling<sup>2</sup>, E. Mauerhofer<sup>1</sup>, M. Willenbockel<sup>2</sup>

<sup>1</sup>Institute of Energy and Climate Research - Nuclear Waste Management and Reactor Safety, Forschungszentrum Jülich GmbH, 52425 Jülich, Germany <sup>2</sup>Central Institute for Engineering, Electronics and Analytics – Electronic Systems, Forschungszentrum Jülich GmbH, 52425 Jülich, Germany

# Introduction

For non destructive characterization of nuclear waste detailed information about massive and dense structural components are needed from radiography to improve analytical results.

Setup

#### **Detector Design**

- Commercial X-Ray detector (PerkinElmer)
- Active area: 40 x 40 cm<sup>2</sup>
- Segmentation: 1024 x 1024 pixels




Fig. 1: Opened detector without scintillator. 1 Active area made from aSi. 2 Electronics

### Scintillator

- General purpose plastic scintillator
- EJ-260 Eljen Technology
- Thickness: 3 mm



Fig. 2: General purpose plastic scintillator used for fast neutron detection via recoil protons.

# **Neutron generator**

- Commercial generator Genie16GT (Sodern)
- D-T fusion for 14 MeV neutrons
- Flux determination with monitoring foils (Al, Au)
- Distance source to foils: 30 cm
- Activity measurement with HPGe detector

• Fast neutron source strength:

Reaction	A in Bq	Φ in cm <sup>-2</sup> s <sup>-1</sup>	Q in s <sup>-1</sup>
<sup>197</sup> Au(n,γ) <sup>198</sup> Au	2.48	28.1·10¹	3.2.106
<sup>197</sup> Au(n,2n) <sup>196</sup> Au	0.10	12.2·10 <sup>3</sup>	1.4·10 <sup>8</sup>
$^{27}$ Al(n, $\alpha$ ) $^{24}$ Na	0.38	15.5·10 <sup>3</sup>	1.7-108

 $Q_{fast} = 15.7 \pm 2.6 \cdot 10^7 \,\text{n/s}$ 

**Tab. 1**: Analysed reactions with corresponding activities, neutron fluxes and source strengths.

## **Experiment**

- Neutron generator within 10-20 cm PE shield
- Distance source to detector: 42 cm
- Samples on lift table




Fig. 3: Experimental setup with samples of PE and lead. Shielding lid is removed.

# **Image Analysis**

### Setup

- Pb brick and PE cylinderDistance source detector: 42 cm
- Average of 450 frames, each 2 s





Fig. 4: Photograph of the setup for the following radiographs.

# Smooth

- Profiles of areas with and without objectsSet outliners to the average of the surrounded
- pixels



**Fig. 5**: Profiles of the area without object (red) and including the object (blue).

# Profile CorrectionSet area without objects as background







Fig. 6: Profiles of the smoothed radiograph.

# Signal AnalysisFit Gaussian distribution to histogram from

region of interest








**Fig. 7**: a) Profiles of the corrected radiograph. b) histogram of the region of interest from areas seen in c). d) Mean values of the Gaussian fits.

# Results

# Calibration

- Radiographs of well known test samples
   Size: 5 x 8 x 10 cm<sup>3</sup>
- Size: 5 x 8 x 10 cm<sup>3</sup>
  Al, C, Fe, Pb, W, concrete, PE
- Al, C, Fe, PbPE as reference
- Combination of two samplesAnalysis as shown before




**Fig. 9**: a) Radiograph of Pb (red) and PE (black) sample. b) correlation between measured and calculated signal attenuation.

# Summary

First radiography with test samples successful, despite ow detector efficiency and neutron intensity

Discrimination between light and heavy objects

Correlation between detector signal and absorption properties







Fig. 10: Photograph of the experimental setup and corresponding radiograph of an eye bolt M52.


y iy 0 18=2 ≥ 3 00 4 2 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 2 ● 0 21 z

# **New Scintillator**

Fig. 8: Test samples (AI, Graphite, Fe, Pb, W, concrete, PE).


- Stack of scintillating fibres for increased neutron
- conversion efficiencyType: SCSF-3HF(1500)MJ from Kuraray
- Thickness: 10 mm, diameter: 1 mm





**Wavelength Shifting Fibres Detector** 

- Prototype detector
- Plastic with ZnS as scintillator
   X X crossed WLS fibros
- X-Y crossed WLS fibresFibre readout with PMT
- TDC coincidences for position reconstruction
   Active area: 4 x 4 cm² (16 x 16 fibres)








Fig. 12: a) Photomultiplier used for the detector. b) Close-up from the crossed fibres. c) Detector in light tight housing.

**Outlook**