

EVANESCENT WAVE DYNAMIC LIGHT SCATTERING AT A LIQUID-LIQUID-INTERFACE

Internship Report: from data evaluation to a new measurement setup

20. APRIL 2021 | SAMUEL MONTER

EVANESCENT WAVE DYNAMIC LIGHT SCATTERING

Measurement Setup



- Particles scatter light from evanescent wave

$$I = I_0 e^{-kz} \text{ with } \frac{\kappa}{2} = \frac{2\pi}{\lambda} \sqrt{(n_1 \sin(\alpha_i))^2 - n_2^2}$$

- 3 axis goniometer enables resolution of diffusion normal and parallel to the interface

$$Q_{\parallel} = \frac{2\pi n_2}{\lambda} \sqrt{1 + \cos(\alpha_r^2) - 2 \cos(\alpha_r) \cos(\theta)}$$

$$Q_{\perp} = \frac{2\pi n_2}{\lambda} \sin(\alpha_r)$$



EVANESCENT WAVE DYNAMIC LIGHT SCATTERING

Data Evaluation

- Measure: mixed homodyne and heterodyne detection

$$g_2 = 1 + 2C_1 g_1(t) + (C_2 g_1(t))^2$$
$$\text{with } C_1 = C_2 - C_2^2 \quad C_2 = 1 - \sqrt{1 - A}$$
$$g_1 = (1 - B_1) \exp(-\Gamma t) + B_1$$

- Extract: Γ from a fit to g_2
 - Model for short times for Γ


$$\Gamma = Q_{\parallel}^2 \langle D_{\parallel} \rangle (\kappa) + \left(Q_{\perp}^2 + \frac{\kappa^2}{4} \right) \langle D_{\perp} \rangle (\kappa)$$

- Decay is single exponential for short times

EVANESCENT WAVE DYNAMIC LIGHT SCATTERING

Data Evaluation – so far


- Choose start and end point for fit
 - Perform fit $\rightarrow \Gamma$
 - Cut off data points for long times
- end point vs. Γ plot
- Choose a Γ from constant region



EVANESCENT WAVE DYNAMIC LIGHT SCATTERING

Data Evaluation – Extension


- Choose start and end point for fit
 - Perform fit $\rightarrow \Gamma$
 - Cut off data points for long times
 - Cut off data points for long times
- End point vs start point vs Γ contour plot
- Choose a Γ from constant region
 - Benefits:
 - Less human bias
 - Automated stack execution




EVANESCENT WAVE DYNAMIC LIGHT SCATTERING

Measurements on hollow sphere

- Determine $\langle D_{\parallel} \rangle$ from the slope of a linear fit of $Q_{\parallel}^2$ vs Γ




$$\Gamma = Q_{\parallel}^2 \langle D_{\parallel} \rangle (\kappa) + \left(Q_{\perp}^2 + \frac{\kappa^2}{4} \right) \langle D_{\perp} \rangle (\kappa)$$



NEW SETUP: LIQUID-LIQUID-INTERFACES


New challenges

- Choice of the system
 - Immiscible liquids
 - $n_1 < n_2$
 - Stability against irreversible adsorption at the interface
 - Negatively charged PTFE particles in H_2O + toluene
- Setup alignment
 - Alignment of empty cell
 - Precise filling of the cell




NEW SETUP: LIQUID-LIQUID-INTERFACES

Cell Alignment




NEW SETUP: LIQUID-LIQUID-INTERFACES

Cell Alignment – x,y,z




NEW SETUP: LIQUID-LIQUID-INTERFACES

Cell Alignment – x,y,z




NEW SETUP: LIQUID-LIQUID-INTERFACES

Cell Alignment – x,y,z





NEW SETUP: LIQUID-LIQUID-INTERFACES

Cell Alignment – x,y,z



NEW SETUP: LIQUID-LIQUID-INTERFACES


Cell Alignment – α_i, α_r



$$\Delta\alpha = \frac{1}{2}(\alpha_i - \alpha_r)$$


NEW SETUP: LIQUID-LIQUID-INTERFACES

Cell Alignment - h





NEW SETUP: LIQUID-LIQUID-INTERFACES

Cell Alignment





NEW SETUP: LIQUID-LIQUID-INTERFACES

First results




NEW SETUP: LIQUID-LIQUID-INTERFACES

First results



NEW SETUP: LIQUID-LIQUID-INTERFACES

First results



NEW SETUP: LIQUID-LIQUID-INTERFACES

Outlook


- Necessary improvements:
 - Cell with higher optical quality
 - Larger cell
- New available experiments:
 - EWDLS at liquid-liquid interfaces
 - Kinetics of nanoparticle interface adsorption
 - Diffusion of particles adsorbed in the interface

THANKS FOR YOUR KIND ATTENTION!

Lets have a nice discussion!


EVANESCENT DYNAMIC LIGHT SCATTERING


Measurements on hollow sphere



MICROFLUIDIC SETUP


1. Vacuum supply
2. Power supply
3. USB
4. Pressure supply
5. Line to the sample cell





NEW SETUP: LIQUID-LIQUID-INTERFACES

Measurement



NEW SETUP: LIQUID-LIQUID-INTERFACES

First results

